

PRODUCT CATALOG for AskDerm.com
http://www.askderm.com
Generated on 2/8/2011

Age Sciences Product Line

Age Sciences - Personal Microderm Tool Kit

Product ID: 39151



Achieve amazing anti-aging results at home with Personal Microderm System. It works just like a professional microdermabrasion treatment without the hassle or cost. The hand-held device uses the same aluminum oxide crystals as the ten thousand dollar machines found at your doctor's office. It smooths away dead skin cells, diminishing age spots and other signs of aging.

The system includes: 1 personal microderm, 6 exfoliating discs (2 green and 4 blue), 1 electric cord (Standard U.S. Outlet), 1 facial cap, 1 body cap, 1 reusable filter, Instructional CD and Pamphlet

All Sales Final - No Returns or Exchanges. If product is defective at the time of receipt call 1-800-393-3848.

Retail Price: \$179.99

Usage:

Step 1: Assembly

Place disc on motor shaft. Simply push on tightly. To remove; Pull Off. Important: There are 2 disc sizes. One for the body and one for the face. See recommended usage below. Place the cap over the disc. Screw 1/4 turn and seal tight. There must be a seal in order to create the vacuum and suction. Connect cord to unit and power supply and plug in.

STEP 2: Wash Your Skin

Cleanse and Wash. Remove all oil, makeup, and dirt. Make sure your hair is out of the way and covered or pulled back off of your face.

STEP 3: Technique

Pull skin tight for a smooth surface. Place your fingers on the skin and stretch the skin so that it is tight where the tip will make contact with the skin. Keep the unit moving constantly over your skin. Do not hover in one spot.

Be sure the cap is flush with your skin, maintaining even contact. This will ensure a proper suction. The motor sound will change when suction occurs.

Small Disc: Use the small disc around tight spots on the nose or where your bones will not allow suction. The large disc should only be used on the body and not on the face.

STEP4: Treatment

Use light pressure with early treatments. As you become skilled with this tool you can become more aggressive and use more pressure. Be careful not to grind to deep. Move the disc in upward and horizontal motions until you have treated all areas of the surface. If you have trouble spots take a little more time on these areas. Be careful not to press to hard or stay to long in one place. Keep the disc moving over the surface or your skin.

Warning Warning Warning Warning Warning Warning Warning Warning Warning

Hovering in one place too long may cause irritation and result in a scab. Treating trouble spots. Both time and pressure in trouble spots will determine how quickly you see results. Blue Disc is Moderate. Green Disc is Coarse for more aggressive treatments.

STEP 5: Cell Renewal

Use our recommended products or your own. The Cell rejuvenation and renewal serums and creams are essential to the healing process and recovery. Use them daily in between treatments.

- Apply Treatment Toner. This neutralizes the pH balance in your skin and prepares the skin for the Recovery Moisturizer. Wait 5 minutes.
- Apply a recovery Moisturizer or Cell Renewal Lotions or serum.
- Use sunscreen for at least 3 days after your treatment; your skin will be extra sensitive.
- Cosmetics may be applied after the treatment. Apply these each day until your next treatment. One week.

STEP 6: Follow-Up and Maintenance

Clean the white filter and residue out of the cap and motor shaft. Clean with a damp cloth in and around motor shaft.

Replace discs after every two to three treatments. For maximum effectiveness, the abrasive on the discs must be fresh.

Recommended usage: Once a week. Note: The Green Disc has a Coarse Grit for more aggressive treatments. The Blue has Moderate Grit. Due to the personal nature of this tool we accept no returns. 90 day warranty on parts and motor.

Age Sciences - Replacement Disc Blue Large

Product ID: 39387

Replacement Disc for Personal Microderm Device.



Ideal for those with age spots, skin discolorations, and other imperfections who want smooth, even-toned skin.

All Sales Final - No Returns or Exchanges.

Retail Price: \$13.00

Usage:

Use on your face, elbows, hands, feet, or wherever you have age spots or fine lines and wrinkles. Works well on surface scars.



Age Sciences - Replacement Discs Blue Mix

Product ID: 39640

Replacement Disc for Personal Microderm Device.

Ideal for those with age spots, skin discolorations, and other imperfections who want smooth, even-toned skin.

Personal Microderm Blue Replacement Discs 6 pack small tips. Sensitive Grit aluminum oxidized coating on tip.

Retail Price: \$13.00

Usage:

Use on your face, elbows, hands, feet, or wherever you have age spots or fine lines and wrinkles. Works well on surface scars.



Age Sciences - Replacement Discs Blue Small

Product ID: 39393

Personal Microderm Blue Replacement Discs 6 pack small tips. Sensitive Grit aluminum oxidized coating on tip.

All Sales Final - No Returns or Exchanges.

Retail Price: \$13.00

Usage:

Use on your face, elbows, hands, feet, or wherever you have age spots or fine lines and wrinkles. Works well on surface scars.



Age Sciences - Replacement Discs Green Large

Product ID: 39394

Age Sciences Green Replacement Disc 6 pack, large. Moderate Grit aluminum oxidized coating on tip.

All Sales Final - No Returns or Exchanges.

Retail Price: \$13.00

Usage:

Use on your face, elbows, hands, feet, or wherever you have age spots or fine lines and wrinkles. Works well on surface scars.

Age Sciences - Replacement Discs Green Mix

Product ID: 39641

Age Sciences Green Replacement Disc 3 pack, large. Moderate Grit aluminum oxidized coating on tip.

Age Sciences Green Replacement Disc 3 pack small tips. Moderate Grit aluminum oxidized coating on tip. Ideal for those with age spots, skin discolorations, and other imperfections who want smooth,



even-toned skin.

Retail Price: \$13.00

Usage:

Use on your face, elbows, hands, feet, or wherever you have age spots or fine lines and wrinkles. Works well on surface scars.



Age Sciences - Replacement Discs Green Small

Product ID: 39395

Age Sciences Green Replacement Disc 6 pack small tips. Moderate Grit aluminum oxidized coating on tip. Ideal for those with age spots, skin discolorations, and other imperfections who want smooth, even-toned skin.

All Sales Final - No Returns or Exchanges.

Retail Price: \$13.00

Usage:

Use on your face, elbows, hands, feet, or wherever you have age spots or fine lines and wrinkles. Works well on surface scars.



Age Sciences - Replacement Discs Red Large

Product ID: 39396

Personal Microderm Red Replacement Discs 6 pack of large sized tips. Coarse Grit aluminum oxide coated tip.

All Sales Final - No Returns or Exchanges.

Retail Price: \$13.00

Our Price: \$10.00

Usage:

Use on your face, elbows, hands, feet, or wherever you have age spots or fine lines and wrinkles. Works well on surface scars.



Age Sciences - Replacement Discs Red Small

Product ID: 39397

Personal Microderm Red Replacement Discs 6 pack, small. Coarse Grit aluminum oxide coated tip.

All Sales Final - No Returns or Exchanges.

Retail Price: \$13.00

Usage:

Use on your face, elbows, hands, feet, or wherever you have age spots or fine lines and wrinkles. Works well on surface scars.